

Model J-4002 Stereo iPad Connected to AUX Input

Bose Acoustic Wave Music System II
AM/FM/CD Player

Speakers Built-In to Be Heard in Area Near Console

Optional iPad (not included with stereo) Connected to Wi-Fi Network Playing Continuous Music

Pandora
Spotify
iTunes
Apple Radio
MP3 Music

Music from Bose Stereo

Aux Input

3.5mm Stereo Plug

Mixer / Amplifier Adjusts Sound for Patient

Microphone Positioned at Operator Console Enables Tech to Speak to Patient

Shielded Scan Room

Optional RF Filter Inserts

RF Penetration Panel

MRI Magnet

Optional Noise Reduction Headphone

Headphones Hung On Plastic Clips When Not In Use

Patient Acoustic Interface

